

### DNA Extraction Kit

Cat. #BT06AB

Store at 15-25°C

### Notes before starting

- All centrifugation steps are carried out at room temperature (15–25°C)
- Heat a water bath or heating block to 65°C  
(for animal tissues heat another one to 70°C)

### Genomic DNA Extraction Protocol

#### 1. Sample preparation

##### a) Cultured Cells

Harvesting monolayer adherent cells after cell count (do not use more than  $5 \times 10^6$  cells)

**Aspirate** the culture medium

**Wash** cells using PBS (1X) then aspirate again

**Detach** cells using trypsin (0.10 – 0.25%) or a cell scraper

**Centrifuge** for 5 min at 300 x g. Discard the supernatant and avoid disturbing the cell pellet

**Resuspend** cells in 200 µl PBS (1X)

##### b) Whole Blood

**Collect** the whole blood in the presence of anticoagulants such as EDTA (10 µl of 0.5M EDTA/ml Blood)

##### c) Animal Tissue

**Cut** ≤25 mg of animal tissue sample and transfer it to a 1.5 ml microcentrifuge tube.

**Homogenize** using mechanical homogenization

2. Add 20 µl Proteinase K solution into a 1.5 ml microcentrifuge tube.  
**Optional:** If RNA interferes with your downstream applications, add 4 µl RNase A (100 mg/ml), mix by vortexing, and incubate for 2 min at room temperature
  3. Transfer 200 µl of the sample to the tube.
  4. For animal tissue:
    - a. Incubate at 65°C for 15 minutes. Mix by inversion occasionally until the cells are completely lysed.
    - b. Add 200 µl of Lysis Buffer to the tube. Mix thoroughly by vortexing or pipetting.
    - c. Incubate at 70°C for 10 minutes.
- For other samples:
- a. Add 200 µl of Lysis Buffer to the tube. Mix thoroughly by vortexing or pipetting to obtain a uniform suspension (A homogeneous mixture is essential for efficient lysis).  
**Note:** Solution may become viscous after addition of Lysis buffer
  - b. Incubate at 65°C for 15 minutes. Mix by inversion occasionally until the cells are completely lysed.
5. Add 200 µl of absolute ethanol (96–100%), Mix thoroughly by pipetting or vortexing for 15 seconds.
  6. Transfer all volume (~620 µl) to the spin column, Centrifuge for 1 min at 10,000 x g. Discard flow through and the collection tube.
  7. Add 500 µl of Wash Buffer I, centrifuge for 1 min at 10,000 x g. Discard the flow through.
  8. Add 500 µl of Wash Buffer II, centrifuge for 3 min at 10,000 x g.
  9. Discard the flow through.  
**Optional:** Repeat steps 8 – 9
  10. Re-spin for 1 min. at max speed (~16,000 x g) to dry the membrane completely
  11. Discard the collection tube and insert the spin column into a recovery tube.
  12. Transfer 25-100 µl of elution buffer or Nuclease-Free Water directly into the spin column membrane.
  13. Incubate the spin column at room temperature for 2-5 min.
  14. Centrifuge for 1 min at 8,000 x g.  
**Optional:** Repeat elution in the same or new tube
  15. Store DNA at -20°C.



For more details,  
Scan me

## Blutruve

9 Kamel Almahrouqi Street

Al Haram, Egypt

+02 237550489

info@blutruve.com